



SHAKING WATER BATH

SHAKING WATER BATH

Microprocessor based, digital shaking baths for gentle heating with agitation for your samples. Easy to clean stainless steel chamber, equipped with clamps that keep samples in place avoiding the risk of spills.

100 SERIES SHAKING WATER BATH



Microprocessor based keypad for temperature/speed controls
 Digital display
 Easy to clean stainless steel internal chamber
 Equipped with clamps to keep samples in place
 It is featured with Spring Bottle rack, Alarm System, Adjustable Thermostat and Electronic Calibration

SPECIFICATIONS

Model	BBSW-101	BBSW-102	BBSW-103
Capacity	50 L		
Maximum Configuration	6 x 500ml or 3-4 x 1000ml		
Rotation Range	0~220 rpm	0~300 rpm	0~220 rpm (orbital) 0~300 rpm(Reciprocating)
Temperature Range	RT+5°C to 99°C		
Time Range	0-999 h		
Platform Dimension	380x300 mm		
Outer Dimension (WxDxH)	705x515x470 mm		
Weight	38 kg		46 kg
Method	30 mm	20 mm	30 mm; 20 mm
Thermo Lid	Orbital	Reciprocal	Orbital and Reciprocal
Power	1600 W		
Power Supply	220V, 50Hz		

200 SERIES SHAKING WATER BATH



Microprocessor based keypad for temperature/speed controls
 Digital display
 Easy to clean stainless steel internal chamber
 Equipped with clamps to keep samples in place
 It is featured with Spring Bottle rack, Alarm System, Adjustable Thermostat and Electronic Calibration

SPECIFICATIONS

Model	BBSW-201	BBSW-202	BBSW-203
Capacity	32 L		
Maximum Configuration	6 x 500ml or 2 x 1000ml		
Rotation Range	0~220rpm	0~300 rpm	0~220rpm (orbital) 0~300 rpm(Reciprocatin)
Temperature Range	RT+5°C to 65°C		
Time Range	0-999 h		
Platform Dimension	430x340 mm		
Outer Dimension (WxDxH)	705x455x440 mm		
Weight	38 kg		46 kg
Method	30 mm	20 mm	30 mm; 20 mm
Thermo Lid	Orbital	Reciprocal	Orbital and Reciprocal
Power	400 W		
Power Supply	220V, 50Hz		



Biolab Scientific Limited

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada
Email: contact@biolabscientific.com | Phone: +1 416 639 2438 | Website: www.biolabscientific.com