



**MICROSCOPIO BIOLOGICO TRINOCULAR
MARCA ADVANCED OPTICAL
REFERENCIA :A12.0701-F
MODELO : BM1000 TR**



A12.0701 Biological Microscope		A	C	E	F
Head	Monocular Head, Inclined at 30°, 360° Rotatable	•	-	-	-
	Seidentopf Binocular Viewing Head, Inclined at 30°, 360° Rotatable, Interpupillary 48-75mm	-	•	•	-
	Seidentopf Trinocular Viewing Head, Inclined at 30°, 360° Rotatable, Interpupillary 48-75mm	-	-	-	•
Eyepiece	WF10x/18	•	•	•	•
Objective	Achromatic objective 4x,10x,40x,100x	•	•	•	•
Nosepiece	Quadruple nosepiece	•	•	•	•
Focusing	Coaxial Coarse and Fine Adjustment, Fine Division 0.002mm, Coarse Stroke 37.7mm per Rotation, Fine Stroke 0.2mm per Rotation, Moving Range 20mm	•	•	•	•
Stage	Double Layers Mechanical Stage132×142mm/75×40mm	•	•	•	•
Condenser	Abbe Condenser NA1.25 with Aperture Diaphragm	•	•	•	•
Illumination	3W LED	•	-	-	•
	Rechargeable 3W LED	-	-	•	-
	12V/20W Halogen Lamp, Brightness Adjustable	-	•	-	-
Optional Accessories		Item No.			
Eyepiece	WF16x/12				
Objective	Plan Achromatic Objective 4x,10x,40x,100x				
	Infinite Plan Achromatic Objective 4x,10x,40x,100x				
Illumination	Kohler Illumination				
Condenser	Dark Field Condensor				
Attachment	Photo Attachment				
	Video Attachment				