



CABINA DE FLUJO LAMINA



LCB-010V Black
Vertical Air Flow Type

ESPECIFICACIONES

Model	LCB-0101V	LCB-0121V	LCB-0151V	LCB-0181V	
Vertical Air Flow Type					
Dimensions (WXDXH mm)	Interior	890 x 550 x 595	1190 x 550 x 595	1490 x 550 x 595	1790 x 550 x 595
	Exterior	1050 x 800 x 1170	1350 x 800 x 1170	1650 x 800 x 1170	1950 x 800 x 1170
Model	LCB-0102H	LCB-0122H	LCB-0152H	LCB-0182H	
Horizontal Air Flow Type					
Dimensions (WXDXH mm)	Interior	890 x 480 x 650	1190 x 480 x 650	1490 x 550 x 595	1790 x 550 x 595
	Exterior	1050 x 800 x 1170	1350 x 800 x 1170	1650 x 800 x 1170	1950 x 800 x 1270
Air Volume	980 m ³ / h	1350 m ³ / h	1700 m ³ / h	2030 m ³ / h	
HEPA Filter	Class 100, Purity 0.3 μm Particle Removal 99.99 %				
Pre-filter	3~30 μm Particle Removal, Built-In Back Panel				
Air flow velocity	0.3 ~ 0.6 m/sec				
Controller	9 Stage Velocity Control, Total Using Timer, UV/FL Lamp On/Off				
Material	Interior	Stainless Steel 304			
	Exterior	Epoxy Powder Coated Electro Galvanized Steel			
Noise Level	56 ~ 60 dBA				
Light	UV 20 W × 1 EA, FL 20W × 1 EA		UV 40 W × 1 EA, FL 40W × 1 EA		
Utility Device	Air Cock, Gas Cock, 110 / 220V Consent				
Door	Tempered Safety Glass Sliding Door				
Color (Standard)	Blue & White for Overseas Sales / Black & White for Domestic Sales <small>*can be changed vice versa upon request</small>				
Safety	Over Temp. Over Pressure and Over Current Protector, Safety Release Valve)				



ADVANCED INSTRUMENTS S.A.S

La ingeniería puesta al servicio de su Empresa.

Safety	Over Temp. Over Pressure and Over Current Protector, Safety Release Valve)			
Voltage	110 V, 60 Hz or 220V,50/60Hz, 1 Phase			
Power Consumption (220V)	3 A	3 A	3 A	6 A

*Specifications can be changed without prior notice for quality improvement
*Custom made order is available

Estimated Packing Details				
<i>*Actual packing details may be different. This is for your reference only</i>				
Model	LCB-0101V	LCB-0121V	LCB-0151V	LCB-0181V
	LCB-0102H	LCB-0122H	LCB-0152H	LCB-0182H
Dimensions (WXDXH mm) (Approx, including, carton)	-	-	-	-
Net Weight Approx	-	-	-	-
Gross Weight Approx	-	-	-	-

Optional Accessories / Options				
Model	LCB-0101V	LCB-0121V	LCB-0151V	LCB-0181V
	LCB-0102H	LCB-0122H	LCB-0152H	LCB-0182H
Movable Stand	LST-0100	LST-0120	LST-0150	LST-0180
ULPA Filter	o	o	o	o
IQ/OQ Validation Documentation at Factory	o	o	o	o