



CONGELADOR



LBF-4015U

MODEL		LBF-4010U	LBF-4015U	LBF-4020U	LBF-4025U
Type		Upright			
Dimensions (W x D x H mm)	Interior	600 x 665 x 960	600 x 665 x 1280	750 x 665 x 1280	750 x 665 x 1280
	Exterior	835 x 855 x 1600	835 x 855 x 1925	985 x 855 x 1925	985 x 975 x 1925
Capacity		369 ℓ	492 ℓ	615 ℓ	730 ℓ
Compressor		1/2 HP	3/4 HP	1 HP	3/2 HP
Insulation		High Density Polyurethane (Thickness : 80 mm)			
Temperature	Range	10°C ~ -35°C			
Controller		Microprocessor Digital Controller			
Display		4 Digit LED Display & Signal Lamps			
Material	Interior	Stainless Steel Polished			
	Exterior	Epoxy Powder Coated Steel			
	Inner Door	Inner Plastic Insulated Door			
Door		Door Key Lock, Magnetic Packing			



ADVANCED INSTRUMENTS S.A.S

La ingeniería puesta al servicio de su Empresa.

Shelves	2 EA	3 EA		
Safety	High/Low Temp. Alarm, Electric Current Leakage Breaker, Audible and Visible Alarm			
Voltage	220/230 V 50/60 Hz, 1 Phase (Option: 110V, 60Hz)			
Power Consumption (220V)	10 A	10 A	12 A	12 A

*Specifications can be changed without prior notice for quality improvement
*Custom made order is available

Packing Dimensions Approx. including carton (W×D×Hmm)	x x	x x	x x	x x
Net Weight Approx				
Gross Weight Approx				