



CONGELADORES DE LABORATORIO



LLF-401SR



LLF-402SR



LLF-403SR-H



LLF-403SR



LLF-404SR



LLF-404SR-H



ADVANCED INSTRUMENTS S.A.S

La ingeniería puesta al servicio de su Empresa.

MODEL		LLF-401SR	LLF-402SR	LLF-403SR	LLF-404SR
Door	Solid Swing Door	1 EA			2 EA
MODEL				LLF-403SR-H	LLF-404SR-H
Door	Solid Half Door				4 EA
Type		Bench Top	Floor Standing		
Dimensions (W x D x H mm)	Interior	460 x 500 x 630	560 x 660 x 800	620 x 700 x 1410	1140 x 700 x 1410
	Exterior	560 x 606 x 850	683 x 775 x 1490	740 x 805 x 2000	1260 x 805 x 2000
Capacity		105 ℓ	295 ℓ	612 ℓ	1125 ℓ
Compressor		1 / 8 HP	1 / 4 HP	1 / 4 HP	1 / 2 HP
Insulation		High Density Polyurethane (Thickness : 60mm)			
Temperature	Range	0°C ~ -25°C			
	Uniformity	±3°C at -20°C			
Controller		Microprocessor Digital Controller			
Display		4 Digit LED Display & Signal Lamps			
Material	Interior	Stainless Steel 304			
	Exterior	Epoxy Powder Coated Steel			
Air Circulation		Turbo Fan with Automatic Stop/Run during Door Open			
Shelves		2 EA	3 EA	6 EA	12 EA
Safety		High/Low Temp. Alarm, Electric Current Leakage Breaker, Audible and Visible Alarm			
Voltage		220V, 50/60Hz, 1 Phase (Option: 110V, 60Hz)			
Power Consumption (220V)		4 A			