



ADVANCED INSTRUMENTS S.A.S.

La Ingeniería puesta al servicio de su empresa
NIT: 830.101.830-1 IVA REGIMEN COMUN

IMPORTADORA Y COMERCIALIZADORA
DE EQUIPOS PARA CONTROL DE CALIDAD
ENSEÑANZA Y ANALITICOS

DESTILADOR DE AGUA BÁSICO S / S

MODELO LWD

MARCA LABTECH



LWD-104S Black & STS

LWD-204SF Black & STS



Model		LWD-104S LWD-204SF	LWD-108S LWD-208SF
Dimensions (WXDXH mm)		330 x 290 x 580	330 x 365 x 640
Capacity		3.0~3.5 ℓ / hour max	6.0~7.0 ℓ / hour max
Heater		3 kW	6 kW
Cooling Water		Approx. 30 ~ 50 ℓ	Approx. 60 ~ 80 ℓ
Material	Inner	Stainless Steel	
	Outer	Stainless Steel Hair Lined	
	Heater	Stainless Steel Pipe	
Controller		Hydraulic Thermostat with Main Power Switch	
Sealing		Silicone Gasket	
Operation		Fully Automatic Operation	
Indication Lamp		Distill, Full	
Color (Standard) <i>*can be changed vice versa upon request</i>		Blue & STS for Overseas Sales / Black & STS for Domestic Sales	
Safety		Low Water & Over Temperature Protection	
Voltage		220V 50/60Hz, 1Phase	
Power Consumption (220V)		14 A	28 A

*Specifications can be changed without prior notice for quality improvement

*Custom made order is available

Estimated Packing Details

**Actual packing details may be different. This is for your reference only*

Model		LWD-104S LWD-204SF	LWD-108S LWD-208SF
Dimensions (WXDXH mm)		400 X 380 X 700	450 X 400 X 700
Net Weight	Approx	12	14
Gross Weight	Approx	13	20

Descriptions

Model	Descriptions	
LWD-104S	S/S Basic Water Still	3.0~3.5ℓ/Hour
LWD-204SF	S/S Basic Water Still/ w Prefilter	3.0~3.5ℓ/Hour
LWD-108S	S/S Basic Water Still	6.0~7.0ℓ/Hour
LWD-208SF	S/S Basic Water Still/ w Prefilter	6.0~7.0ℓ/Hour