



ADVANCED INSTRUMENTS S.A.S.

La Ingeniería puesta al servicio de su empresa
NIT: 830.101.830-1 IVA REGIMEN COMUN

IMPORTADORA Y COMERCIALIZADORA
DE EQUIPOS PARA CONTROL DE CALIDAD
ENSEÑANZA Y ANALITICOS

DESTILADOR DE AGUA AUTOMÁTICO DOBLE

MODELO LWD

MARCA LABTECH



LWD-3005D Black & White



LWD-3005D Blue & White

Model	LWD-3005D	LWD-3010D
Dimensions (WXDXH mm)	570 x 600 x600	680 x 720 x650
Capacity	3.0 ~ 3.5 ℓ / hour max	6.0 ~ 7.0 ℓ / hour max
Heater	3 kW x2EA	6 kW x 2 EA
Cooling Water	Approx. 35 ~ 60 ℓ	Approx. 60 ~ 90 ℓ
Material	Interior Stainless Steel Boiler Stainless Steel Storage Tank Exterior Stainless Steel Hair Lined	
Color (Standard) <i>*can be changed vice versa upon request</i>	Blue & White for Overseas Sales / Black & White for Domestic Sales	
Storage Tank	8 ℓ	14 ℓ
Operation	Fully Automatic Operation	
Indication Lamp	Distill, Full	
Safety	Low Water & Over Temperature Protection, Electric Leakage Circuit Breaker	
Voltage	220V 50/60Hz	
Power Consumption (220V)	27 A	55 A

*Specifications can be changed without prior notice for quality improvement

*Custom made order is available

Estimated Packing Details

**Actual packing details may be different. This is for your reference only*

Model	LWD-3005D	LWD-3010D
Dimensions (WXDXH mm)	700 X 800 X 690	870 X 760 X 950
Net Weight Approx	40	47
Gross Weight Approx	52	60