



ADVANCED INSTRUMENTS S.A.S.

La Ingeniería puesta al servicio de su empresa
NIT: 830.101.830-1 IVA REGIMEN COMUN

IMPORTADORA Y COMERCIALIZADORA
DE EQUIPOS PARA CONTROL DE CALIDAD
ENSEÑANZA Y ANALITICOS

ESTEREO BINOCULAR CON ZOOM MODELO: JSZ 5B MARCA ADVANCED OPTICAL



JSZ - 5 is new "A" designed with special advanced Anti-electrostatic-meet ESD requirement, Anti-mold-can be used in hot, humid environments, and Airtight-prevents contamination from dust, oil and water. It is specifically intended for Industry and Routine research workers. it is wide field of view and good image quality. It is providing a bright object field that is rich in details and great working distance that allows a variety of microscope tools, it is ideal for sample identification and sorting.

Features:

- Binocular Optical head Inclination at 45° viewing, 360° rotatable.
- Tube Adjust range : 55mm-75mm
- Standard Magnification : 7.0X ~ 30X (Without optional lens)
- Eyepiece : WF10X-20mm
- Zoom Objective : 0.7X ~ 3X
- Zoom Ratio : 4.3 : 1
- Working Distance : 97mm
- Visual field : $\Phi 28.6 \sim \Phi 6.6$ mm
- Illuminator : Integrated universal power supply from 110V to 240V, 50Hz/60Hz, Power Rating : 15W.
LED illumination for incident and transmitted lights
(Top Illuminator: 1W, Bottom Illuminator: 0.5W).

Model No. JSZ-5B (for option)

- Zoom Objective : 1X ~ 4.5X
- Eyepiece : WF10x/20mm

Optional Accessories:

- Eyepiece: WF15X/14mm, WF20x/11mm
- Auxiliary lens: 0.75x and 2x
- Stand: Plain stand.