

Placa calefactora y agitador

가열 교반기



Model		LMS-1001	LMS-1002	LMS-1003
Dimensions (WxDXH mm)	Plate	180 x 180		
	Overall	200 x 260 x 110		
Temperature	Range	N/A	60°C to 380°C (Max.)	
Speed	Range	60 to 1500 rpm	N/A	60 to 1500 rpm
Controller		Electronic Solid State Controller		
Material	Top Plate	Ceramic Coated Stainless Steel		
	Body	Anti-Corrosive Powder Coated Steel		
Type		Magnetic Stirrer	Hotplate	Hotplate & Stirrer
Voltage		110V, 60Hz or 220V, 50/60Hz, 1Phase		
Power Consumption (220V)		0.1 A	0.1 A	3 A