



MICROSCOPIO MONOCULAR
REFERENCIA XSP-116NL
MARCA ADVANCED OPTICAL



MONOCULAR MICROSCOPE REFERENCIA EX30/ MODEL XSP-116NL

Specification	Model				
	XSP-116N	XSP-116B	XSP-116F	XSP-116FB	
Viewing Head	Monocular Inclined at 45 degree	V	V	V	V
	360 degree Rotatable		V	V	V
Eyepiece	WF10X	V	V	V	V
Objective	Achromatic 4X, 10X, 40X(s)	V	V	V	V
Stage	Plain stage with slide clips	V	V	V	V
	Stage size: 90mm X 90mm				
Condenser	Single lens with Disc diaphragm	V	V	V	V
	Coarse Focusing Adjustment	V		V	V
Focusing	Coaxial coarse & fine focusing adjustment		V		V
	Fine focusing scale value 0.002mm				
Light Source	Incandescent lamp	V	V	V	V
	20V/20W or 110V/20W				
Optional Accessory	Eyepiece: WF16X, WF20X				
	Objective: 20X, 60X(s)				
	LED illumination				