

OIL-FREE-VACUUM PUMP



I Features

- ❑ Basic Vacuum Pump for laboratory applications.
- ❑ Oil-less and low maintenance.
- ❑ Low noise level during operation.
- ❑ Large rubber feet add stability and act as vibration dampers.
- ❑ Easy replacement of moisture trap and filter.
- ❑ Vacuum level adjustable up to 650 mm Hg.
- ❑ High flow rate.
- ❑ convertible to pressure pump (with optional pressure kit).

I Applications

- Vacuum Filtration
- Vacuum Desiccators and Cabinets
- TLC sprayers
- Vacuum Ovens (partial evacuation of larger units)

Model	LVP-500
Dimension	272 x 121 x 164
Vacuum	650 mmHg
200 mmHg	2 sec
300 mmHg	3 sec
400 mmHg	4 sec
500 mmHg	6 sec
600 mmHg	11 sec
650 mmHg	18 sec
Max Flow Rate	28 l/min
Motor Power	1/8 HP
Port Thread	1/8" NPT
Noise Level	60 dB

* Specifications can be changed without prior notice for quality improvement.

I Features

- ❑ Rechargeable 35 g butane fuel tank.
- ❑ 60 min to 90 min use on full charge.
- ❑ Adjustable flame height.
- ❑ 1,300 °C flame suitable for microbiology.
- ❑ Compact size for Clean benches and Fume Hoods.
- ❑ Wide and Heavy base ensures stability during operation.

Model	LBN-386A
Dimension (W x D x H mm)	115 x 75 x 81
Capacity	Butane Fuel 35 g
Working Time	> 60 min
Dry Weight	410 g

* Specifications can be changed without prior notice for quality improvement.