

Multi-Posición Placa calefactora y agitador

멀티 가열 교반기



| Model | | LMS-3003 | LMS-3006 | | |
|--------------------------|-----------|-------------------------------------|----------|-------------------------|------------------|
| Dimensions (WXDXH mm) | Plate | 180 x 180 x 3 Positions | | 180 x 180 x 6 Positions | |
| | Overall | 600 x 260 x 120 | | 600 x 260 x 120 | 1600 x 520 x 220 |
| Position | | 1 x 3 Positions | | 1 x 6 Positions | 2 x 3 Positions |
| Temperature | Range | 60°C to 300°C | | | |
| Heater | | 500 W x 3 EA | | 500 W x 6 EA | |
| Speed | Range | 60 to 1500 rpm | | | |
| Controller | | Electronic Solid State Controller | | | |
| Material | Top Plate | Ceramic Coated Stainless Steel | | | |
| | Body | Anti-Corrosive Powder Coated Steel | | | |
| Voltage | | 110V, 60Hz or 220V, 50/60Hz, 1Phase | | | |
| Power Consumption (220V) | | 8 A | | 15 A | |