



REFRIGERADOR



LCV-201GR

LCV-202GR

LCV-212GR



LCV-203GR

LCV-204GR



MODEL		LCV-201GR	LCV-202GR	LCV-212GR	LCV-203GR	LCV-204GR
Type		Bench Top	Floor Standing			
Dimensions (W x D x H mm)	Interior	460 x 420 x 630	560 x 660 x 800	500 x 640 x 1400	620 x 700 x 1410	1140 x 700 x 1410
	Exterior	560 x 627 x 850	683 x 796 x 1490	600 x 750 x 1640	740 x 826 x 2000	1260 x 826 x 2000
Capacity		105 ℓ	295 ℓ	421 ℓ	612 ℓ	1125 ℓ
Compressor		1 / 8 HP	1 / 4 HP	1 / 4 HP	1 / 4 HP	1 / 2 HP
Temperature	Range	2°C ~ 10°C				
	Uniformity	± 1°C at 4°C				
Display		4 Digit LED Display & LED Signal Lamps				
Controller		Microprocessor Digital Controller with Battery Back Up System				
Material	Interior	Stainless Steel 304				
	Exterior	Epoxy Powder Coated Steel				
Insulation		High Density Polyurethan (Thickness : 60 mm)				
Air Circulation		Turbo Fan with Automatic During Door Open				
Drawer Shelves		2EA & 1Basket	3 EA	6 EA	6 EA	12 EA
Door		Vacuum Sealed Pair Glass Swing Door with Invisible Heat Film Embedded, Key locker, Adjustable Tension Spring, Magnetic Packing Door				

Safety		High/Low Temp. Alarm, Electric Current Leakage Breaker, Audible and Visible Alarm				
Voltage		220V, 50/60Hz, 1 Phase (Option: 110V, 60Hz)				
Optional Size Available for LCV-204GR	Interior	1370 x 700 x 1440				
	Exterior	1490 x 850 x 2000				
	Capacity	1352 ℓ				
Power Consumption (220V)		2 A	4 A			

*Specifications can be changed without prior notice for quality improvement
*Custom made order is available

Packing Dimensions Approx. including carton (WxDxHmm)	x x	850 x 720 x 1860	900 x 800 x 1640	810 x 880 x 2040	1370 x 960 x 2050
Net Weight	Approx	172	128	150	270
Gross Weight	Approx	187	143	165	285



Optional Accessories	LCV-201GR	LCV-202GR	LCV-212GR	LCV-203GR	LCV-204GR
Thermal Recorder	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Mobile Alarm	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Dual Evaporators	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Dual Compressors	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
7 Days Chart Recorder		<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Stainless Steel Drawer	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Security Box with Lock	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
PC Communication (RS-232C or RS-485)	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
IQ/OQ Validation Documentation at Factory	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>