



**CE**  
Certified

- ❖ Compact Personal Microcentrifuge
- ❖ Maximum RCF up to 2,200 x g
- ❖ Durable Polycarbonate Lid and ABS Resin Body

## Compact and Slim

LabTech Microcentrifuges are compact and economy for personal use

경제적인 가격의 개인용 microcentrifuge

## Choice from Three Models

Choice from three models that spin either 1.5ml or 2.0ml microtubes, 6 adapter for 0.5ml microtube, 8-well tube strip and 0.2ml tube adapter available

다양한 튜브를 사용할 수 있는 모델

## Convenient Operation

Power switch on right side of the body simply close lid to start and open to stop spin

커버를 닫으면 작동되고 열면 중지하는 편리성

## Durable and Stable

New hinged design and latch provide durable and stable operation

내구성이 강한 Polycarbonate 재질의 뚜껑과 ABS 수지 재질 사용

Model	GMC-060	GMC-160	GMC-260
Rotor Type	2.0 ml Tube	0.2 ml Tube	Both
Max Speed	6600 rpm	6300 rpm	Both
Max RCF	2,000 x g	2,200 x g	Both
Material	Lid	Polycarbonate	
	Body	ABS Resin	
	Rotor	ABS Resin	
	Rubber	NBR	
Dimensions	126 (W) x 138 (D) x 100 (H) mm		
Electric Supply	90 ~ 240 VAC Free Voltage, 50/60Hz		

*Model GMC-260 supplied with 6 adapters for 0.5ml microtube, 2.0ml basic rotor for 8 well tube strip*

☒ Custom made order is available  
Specifications can be changed without prior notice for quality improvement.